



Catalog

Electrochemistry accessories.



Your one-stop-shop for accessories and add-ons.

As one of the world's first designers and manufacturers of high-performance electrochemical measurement instruments, BioLogic has forged its place in the international market.

Our comprehensive product portfolio covers cutting-edge scientific products in such diverse applications as electrochemistry, battery testing and fuel cell/material testing.

And supporting these high-precision measurement tools is an extensive line of product accessories, ranging from sophisticated quartz crystal microbalances to electrodes.

So whatever your field of research, we can provide you with high-quality, hard-wearing equipment.

BioLogic accessories: a one-stop-shop to suit your every need.



BioLogic Add-on Instruments.

Rotating electrodes: BluRev	RDE/RRDE systems Tips Cells Enclosure	04-07
------------------------------------	--	--------------

Quartz Crystal Microbalance: BluQCM	Systems Cells Sensors	08-11
--	-----------------------------	--------------

Accessories.

Battery accessories	Coin cell holders Cylindrical cell holders Pouch cell holders Prismatic & pouch cell holders Current collectors Sense adapter module (SAM-50) Redox Flow Batteries (RFB)	12-17
----------------------------	--	--------------

Analytical cells	Small volume cells Large volume cells Multi purpose cells	18-21
-------------------------	---	--------------

Corrosion cells	Standard corrosion cells Avesta cell Flat cells Galvanic cells Coating cells Plate material evaluating cell	22-25
------------------------	--	--------------

Electrodes	Working electrodes Counter electrodes Small reference electrodes Maintenance of reference electrodes Hydrogen reference electrodes	26-31
-------------------	--	--------------

Connection accessories	High temperature cables Glove box cables Multi-electrode investigation cables Connectors Connection kits External device connection Test boxes Faraday cages	32-34
-------------------------------	---	--------------

Spectroelectrochemistry	Spectrometer system Static cell	35-37
--------------------------------	------------------------------------	--------------

Scanning Products Accessories	Probes Video microscope system Cells M470 glove box cables USB-PIO	38-41
--------------------------------------	--	--------------

Material Testing Accessories	Temperature control units Sample holders	42-43
-------------------------------------	---	--------------

BluRev Rotating disk/rotating ring disk electrodes

Steady state efficiency



Exploit the power of EC-Lab with BluRev rotating disk solutions

BluRev is a range of **robust, versatile**, rotating disk (RDE) and rotating ring-disk (RRDE) electrodes ideal for use with BioLogic potentiostats. A range of quick-fit exchangeable electrode tips (multiple materials/diameters) can be easily added to the body.

Driven by EC-Lab software, an industry standard for potentiostat control software, BluRev instruments can be used independently or as part of an integrated setup, where the end-user controls each element of his experiment remotely with total ease.

BluRev and Enclosure: A fully integrated solution

The BluRev's custom designed enclosure has been specifically designed to facilitate the operation



and set-up of the BluRev RRDE. This makes it easy to hold the RRDE body as well as the BioLogic potentiostat cables. The result is a fully integrated, purpose-built rotating electrode system.

Which one do I need? An RDE or RRDE? [What are the advantages of the extra "R" ? Click here](#)

Rotating Electrodes.

BluRev Systems

These modular research instruments demonstrate excellent levels of accuracy, particularly at high rotation rates, and are the only devices on the market to display both target and actual rates of rotation – data that can greatly assist users during experiments by enabling them to monitor rotation rates in real-time.



The RC-10k control unit offers an accuracy of ± 1 rpm over the whole rotational range for precise and fully Reproducible experimental conditions.

The speed can be set manually or remotely by using the analog output of a BioLogic instrument. Any other device with an analog output can also be used.

Specifications

With EC-Lab, you can easily program, perform and analyse ring-disk, Levich and Koutecky-Levich experiments and also use the new EIS element Winf that will directly give you the diffusion coefficient.


See Application note #66 to "EIS measurements on a Rotating Disk Electrode (RDE) Part I: Determination of a diffusion coefficient using the new element Winf"



094-RDE without tip

Rotational range/rpm	100 - 10,000 (9000 for RRDE)
Accuracy/rpm	1 typical over the whole rotational range
Precision/rpm	3 typical over the whole rotational range
Runout/mm	<0.1
Setting resolution/rpm	10 (Manual control) or 25 (Remote control)
Materials of RRDE/RDE	PEEK, Al alloy
Inert gas inlet for shaft corrosion protection/mm	Ø 2
Operating temperature/°C	10 - 40
Power	24 Vdc, 1 A max
Max consumption/W	24
Dimensions/mm	RC-10k: 95x227x178 (HxWxD), RDE/RRDE: 233.6 (length with tip)
Weight/kg	RC-10k: 1.00 RDE/RRDE: 0.36 (without tip)

BluRev RRDE and RDE	Catalog n°	094-RC/RDE	094-RC/RRDE
Content			
RC-10k Rotation controller	094-RC	1	1
Rotating disk electrode (motor, shaft, electrode body with 1 pairs of Ag/C brushes)	094-RDE	1	-
Rotating ring-disk electrode (motor, shaft, electrode body with 2 pairs of Ag/C brushes)	094-RRDE	-	1
DB9 to BNC connector for external control of RC-10k	092-22/1	1	1
1 m BNC/BNC cable	COR28100	1	1
Replacement Ag/C brush	094-RDE-BRUSH	2	4
1 transport case	-	1	1



BluRev Tips

All standard tip bodies (M6 thread) are made of PEEK, but for experiments requiring a high chemical resistance, PTFE versions are also available.

All tips are polished to obtain a final roughness Ra of 50 µm (*Ra = 10 nm for the Boron Doped Diamond tip).


RDE Tips

Disk electrodes for RDE	Catalog n°
Glassy Carbon (3 mm)	094-GC/3
Glassy Carbon (5 mm)	094-GC/5
Glassy Carbon, PTFE body (3 mm)	094-PTFE-GC/3
Glassy Carbon, PTFE body (5 mm)	094-PTFE-GC/5
Silver (3mm), 999 %	094-Ag/3
Aluminum (3 mm), 999 %	094-Al/3
Copper (3 mm), 999 %	094-Cu/3
Nickel (3 mm), 999 %	094-Ni/3
Stainless steel (3 mm), 999 %	094-316L/3
Platinum (2 mm), 999 %	094-Pt/2
Gold (2 mm), 999 %	094-Au/2
Boron Doped Diamond (3 mm) 500 µm diamond layer Boron doping level between 500 and 1000 ppm	094-BDD/3



RRDE Tips

Ring-Disk electrodes for RRDE	Catalog n°
Glassy Carbon Ring & Disk (Nmax = 0.4626)	094-GC-GC
Pt Ring (999 %) - Glassy Carbon Disk (Nmax = 0.4626)	094-Pt-GC



Dimensions/mm: r1=1.5 ; r2 = 2 ; r3 = 3.
N = maximum theoretical collection factor using Albery formula [1].
[1] W. J. Albery and S. Bruckenstein, Trans. Faraday Soc. 62 (1966) 1920.


BluRev Cells

The BluRev RDE/RRDE range is compatible with a wide range of cells – in particular the EL-ELECTRO-80, EL-ELECTRO-80DJ or EL-ELECTRO-150DJ kits.

Please note: if you already possess one of these cells, you will need a special adaptor cap: **094-A-CAP**.

Cells	Catalog n°
PTFE 5 holes cap compatible with BluRev RDE (needed if you already have an EL-ELECTRO cell)	094-A-CAP

BluRev enclosure



BluRev Enclosure	Catalog n°
Protective housing and stand kit for the BluRev	094-ENCL
Contents	
1 stainless steel plate with M6 threaded holes to fix support poles	
2 support poles to hold BluRev RDE and the cell	
(N.B: The clamp and the clamping nut are not included)	
1 clamping flange for the BluRev	
3 half-clamps for the potentiostat cables(all BioLogic cables are supported)	
1 plexiglas protective housing	
Specifications	
Dimensions with protective housing (HxWxD)/mm	287x318x308
Weight (with protective housing)/kg	5.3

BluRev: Background information and theory.

RDE

W_{inf} : direct access to the diffusion coefficient.


Levich and Koutecký-Levich methods¹ are powerful analysis tools used to obtain kinetic electrochemical parameters such as the diffusion coefficient of a redox species in a given medium and the reaction constant. These analyses require potentiodynamic curves at various rotation speeds.

However, fitting impedance measurements made on a redox reaction occurring at a rotating electrode, at only one rotation speed, also enables the direct measurement of the diffusion coefficient.


A PEIS measurement was performed on a 2 mm Pt electrode, using an equimolar solution of $K_3Fe(CN)_6$ and $K_4Fe(CN)_6$ with concentrations of 5 mM in 0.1 M KCl, a BluRev RDE rotating at 2000 rpm.

Using ZFit and the equivalent circuit shown below that contains W_{inf} as a diffusion impedance, we directly obtain the diffusion coefficient of the species of interest, in our case $7.8 \times 10^{-6} \text{ cm}^2/\text{s}$, which is in agreement with the data found in the literature^{2,3}. For more detailed information please see the EC-Lab application note #66⁴.

The considered electrochemical reaction is the following:



The obtained impedance graph is shown above to the right.



RRDE

For RRDE measurements, a bipotentiostat is needed. A bipotentiostat controls the two working electrodes i.e. one channel to control the disk electrode and the other to control the ring. An SP-300 equipped with two channels would be an appropriate instrument for RRDE applications.

Because of the presence of two working electrodes in the same setup, a specific connection mode (to avoid any ground loop trouble) is needed.

This can be achieved in two different ways:

- Isolating the two channels (at least one potentiostat in floating mode),
- Grounding the counter electrode.

The "CE-to-ground" mode of the BioLogic instruments provides this unique capability. Of these two options, the latter is preferred because there is less leakage of current. All of the BioLogic multichannel potentiostats (except cyclers) offer such types of electrode connection. At the disk electrode, the electro-active species are oxidized or reduced according to the applied potential. This new species is detected by reduction or oxidation at the ring, respectively. In a typical experiment, a CV is performed on the disk electrode and a constant voltage is applied on the ring electrode. This is the CV-CA technique that is available in EC-Lab® in the "bipotentiostat" techniques folder.



(1) Application Note #56 "Electrochemical reaction kinetics measurement: the Levich and Koutecký-Levich analysis tools"

(2) A. J. Bard, W. Faulkner, in: *Electrochemical Methods, Fundamentals and Applications*, 2nd Ed., Wiley, New York (2001) 381.

(3) D. R. Lide, H. V. Kehiaian, in: *CRC handbook of thermophysical and thermochemical data*, CRC Press, Inc., Boca Raton, (1994)

(4) Application Note 66 "EIS measurements on a Rotating Disk Electrode (RDE) Part I: Determination of a diffusion coefficient using the new element W_{inf} "

BluQCM Quartz crystal microbalance

Small footprint.
High-sensitivity/reproducibility.



Patented Quick-Lock for **easy setups** and excellent **reproducibility**

The small footprint QCM that packs a punch

The BluQCM QSD is a single channel, compact and modular instrument. Its low footprint and lightweight makes it particularly suitable for crowded labs. It is available as standalone, with temperature control or/and flow control.

The flexibility of the BluQCM QSD is relevant to both electrochemically driven mass weighing applications (electroplating/corrosion, electrode modifications) and the more advanced studies of solid/liquid interfaces (interface rheology, biomolecular interactions).

Quick-lock sensors facilitate set up and increase experimental reproducibility.

Simultaneous measurements can be performed of up to 7 overtones (fundamental included) enabling a complete determination of the layers' properties. Measurements in air, gas, and liquid are possible.



A fast acquisition rate increases mass sensitivity and patented quick-lock sensor cells facilitate set-ups as well as increase the reproducibility of experiments.

All that, and you can control all coupled instruments, including potentiostat/galvanostats, with the same software for eQCM experiments.

BluQCM: A modular approach to QCM

Acoustic sensing principle

The acoustic sensing principle is based on the precise detection of changes on the properties of an acoustic (mechanical) wave traveling through the bulk of the acoustic wave sensor.

quartz crystal induces mechanical oscillations on the quartz due to the piezoelectric effect. A wave is generated and propagated through the sensor and the films attached to it. The resonance frequency of this wave depends on the oscillating mass of the sensor and its adhering layers. When a thin film is attached to the sensor, the properties of the wave change as well, modifying the resonance frequency and amplitude. If the film is thin and rigid, the decrease in frequency is proportional to the mass of the film.

QCM-AWS sensors consist of a thin piece of quartz confined between a pair of metal-based electrodes. An alternating current applied to the

QSD-300: Quartz Crystal



General function	Tracking* and high resolution** at single and multiple overtones
Operation modes	QCM, HFF-QCM, LOVE-SAW
Sensors technologies	See p.3
Cells available	3 - 45 (depending on sensor and cell type)
Liquid volume above sensor/ μL	Optional (QSD-TCU)
Temperature control	Yes
Measurement in air	Up to 7 (up to 13th)
Simultaneous overtones measurements	90x220x260
Dimensions (HxWxD)/mm	3
Weight/ kg	Sensor
Frequency range/MHz	4 - 160
Best frequency resolution/Hz	0.1
Best frequency accuracy/Hz	± 0.5
Max. acquisition rate***/points/s	250
Best mass sensitivity in liquid****/pg/cm ²	50
Best dissipation sensitivity	10^{-7}
Catalog n°	
BluQCM QSD-300	AW-QSD-300

FCU: Flow Control Unit



General function	
Syringe volume/ μL	250 (default)*
Flow rate range for a 250 μL syringe/ $\mu\text{L}/\text{min}$ **	12.5-14500 (Standard) 0.625 - 1062.5 (Smooth)
Dimensions (H x W x D)/mm	195x70x250
Weight/kg	0.75
Catalog n°	
Standard flow control unit	AW-QSD-FCU
Smooth flow control unit	AW-QSD-FCUS

*Other syringe volumes are available upon request, from 12.5 μL to 5000 μL .

** Flow rates depend on the syringe volume. For the standard flow unit, the flow rate change is 0.6250 - 290000 $\mu\text{L}/\text{min}$. For the smooth flow unit, it is 0.0313-21250 $\mu\text{L}/\text{min}$. For more information, please contact your local reseller.

TCU: Temperature Control Unit






General function	
Temperature control range/ $^{\circ}\text{C}$	15 - 45
Temperature stability/ $^{\circ}\text{C}$	± 0.05
Dimensions (H x W x D)/mm	60x220x260
Weight/kg	4.5
Catalog n°	
BluQCM QSD-TCU	AW-QSD-TCU

- [QCM: History and principles](#)
- [QCM: Why measuring at overtones matters](#)
- [QCM: When is the Sauerbrey equation valid?](#)
- [QCM: Measurement principles](#)

QCM: Four articles available on BioLogic Learning Center
Click links to left

BluQCM Cells




eQCM

In-batch eQCM	Flow eQCM	Hermetic Li research in batch
		
AW-BEQ01Q (*AW-BEQATQ for the Air-tight type) *AW-BEQ02Q for 1 inch sensor	AW-FEQ01Q	AW-BEQLIQ

Reference and counter electrodes have to be purchased separately (except for eQCM flow cells where the Pt plate counter electrode is integrated in the lid of the cell).

	Reference electrode		Counter electrode
	Aqueous	Non-aqueous	
In-batch eQCM cells	RE-1B A-012167	RE-7 A-012171	Pt wire 23 cm coiled A-002234
Flow eQCM cells	RE-1S A-012168	RE-7S A-012172	Pt disk integrated in the cell lid

QCM

In-batch	Flow	In-batch probe
		
AW-BQ01Q (14mm sensor) AW-BQ02Q (1" sensor) AW-BQ01HQ (HFF sensor)	AW-FQ01Q (14mm sensor) AW-FQ01HQ (HFF sensor)	AW-PEQ11Q (14 mm sensor)

Sensors

Sensor type	Substrate	Material	Resonant freq./MHz	Finish	Quantity	Catalog N°
14 mm WRAPPED	Cr	Au	10	Polished	10	AW-R10AU10P
	Ti	Au	10	Polished	10	AW-R10AU11P
	-	Cu	5	Polished	10	AW-R5CUP
	-	Al	10	Polished	10	AW-R5ALP
	-	C	10	Polished	10	AW-R10C10P
	-	Pt	10	Polished	10	AW-R10PT10P
	Cr	Au	5	Polished	10	AW-R5AU10P
	Ti	Au	5	Rough	10	AW-R5AU11
	Ti	Au	5	Polished	10	AW-R5AU11P
	Cr	SiO ₂ over Au	5	Polished	10	AW-R5SIO2P


Sensor type	Substrate	Material	Resonant freq./MHz	Finish	Quantity	Catalog N°
1 INCH	Cr	Au	5	Polished	5	AW-R5AU20P
	Ti	Au	5	Polished	5	AW-R5AU21P
	-	Pt	5	Polished	5	AW-R5PT20P
	Cr	Au	9	Rough	5	AW-R9AU20
	Cr	Au	9	Polished	5	AW-R9AU20P
	Ti	Au	9	Rough	5	AW-R9AU21
HFF-QCM	Ti	Au	9	Polished	5	AW-R9AU21P
	Cr	Au	50	-	5	AW-R50AU01H
	Cr	Au	100	-	5	AW-R100AU01H
	Cr	Au	150	-	5	AW-R150AU01H

BluQCM: Background information and theory i

Measurement examples

Electropolymerization of pyrrol

The polypyrrol film was deposited on an Au-coated quartz using cyclic voltammetry (twenty cycles).



The quartz electrode was immersed in an acetonitrile solution (Bu₄NPF₆ 0.2 mol/L) containing a solution of 1 methylpyrrol monomer (0.01 mol/L).

- (1): G. Sauerbrey, Phys. Verh., 1957, 8, 113-114.
- (2): G. Sauerbrey, Z. Phys., 1959, 155, 206-222.
- (3): Application note #13. Section "Apps & literature of EC-Lab division".

Polypyrrol film growth on the quartz working electrode


Fig. 1 represents polypyrrol film growth on the quartz electrode during successive cycles of cyclic voltammetry. The reversibility of the charge transfer in such a polymer film is often dependent on the deposition mode (quasi-reversible in this example). This growth is very regular but tends to slow down during the last cycles. This can be due to an interfacial depletion of the solution in methyl pyrrol species in the layer close to the electrode surface and to a saturation of the working electrode surface area.

QCM measurements during the film growth

Fig. 1 shows the resonant frequency decrease and the resonant resistance increase while the polymer film is growing. Moreover, the variation is dependent on the potential sweep resulting in a pseudo oscillation of frequency and resistance related to successive cycles. This plot can also be made versus potential (see fig. 1).

Fig. 1: overlaid frequency and current vs. E_{we} of the polymer film growth. Scanning at 100mV/s between 0 and 1.018V.

Fig. 2: graphic zoom on one cycle showing the resonant frequency and the current density versus elapsed time⁽³⁾.



Battery Accessories.


Four point cell holders

Four points are better than two

More reliable measurements, higher repeatability levels.

By measuring only the impedance of the cell, a 4-point connection battery holder enables reliable and repeatable measurements. Its design negates measurements of connector and holder-related impedances.

All of our battery holders, with the exception of the CCH model, are built around the 4-point connection design.



Coin Cell Holders



	CCH-1	CCH-8	CCH
Cell max diameter/mm		24	
Cell height/mm	1.6 - 3.2		
Number of channels	1	8	4
Measurement type 4/2 point	4 point	VSP	2 point
		VMP3	MPG-2
		VSP-300	
To be used with	Any instruments	VMP-300	VMP3
		MPG2	
		BCS-805	
		BCS-810	
Climatic chamber compatibility	Yes (-30 to 80°C)	Yes (-30 to 80°C)	No
Catalog n°	096-126	096-128/H Cables must be purchased separately	092-22/14

CCH-8 CABLES



Compatible with	Length	Catalog number
BCS-805	2.5 m long	096-128/C25
BCS-810	5 m long	096-128/C50
VSP / VSP-3e	75 cm long	092-22/24a
VMP-3e		
VSP-300		
VMP-300	2.5 m long	092-22/24b
MPG2		

Cylindrical Cell Holders



BH-1i



CBH-8

2 and 4 mm receptacles are available for the current (power) cables. For voltage (sense) cables, only 2 mm receptacles are available.

4 mm receptacles are available for the current (power) and voltage (sense) cables.



CBH-4

	BH-1i	CBH-4	CBH-8
Cell max diameter/mm	26	60	
Cell min height/mm	0	30	
Cell max height/mm	76	100	
Number of channels	4	4	8
Max current/A	15	32	
Measurement type	4 point	4 point	
Receptacles diameter/mm	2 and 4	4	
To be used with	All instruments		
Max operating T°/°C	60	80	
Size : HxWxD/mm	205x150x95	335x260x150	335x520x150
Weight/kg	0.6	1.9	3.8
Catalog n°	092-22/15	092-C32/4	092-C32/8

More reliable measurements, higher repeatability levels.

Four points are better than two ⁱ

Visit our Learning Center article "[Why four point measurements?](#)" to find out why

Or click [here](#) for the CCH-1 cell holder video.

Pouch Cell Holders



PBH-125

4 POINT measurement



PBH-4

4 POINT measurement

	PBH-125	PBH-150	PBH-4	PBH-8
Min leads separation distance/mm	0			12
Max leads separation distance/mm	110*			44
Number of channels	1	50	4	8
Max current/A	25	50	32	
Measurement type	4 point			
Receptacles diameter/mm	4 (power) 2 (voltage)	6 (power)** 4 (voltage)		4
To be used with	All instruments			
Max operating T°/°C	80	100		80
Size : HxWxD/mm	40x50x210***		135x325x180	135x650x180
Weight/kg	0.2***		1.9	3.8
Catalog n°	092-P25/1	092-P50/1	092-P32/4	092-P32/8

*Measured using the guide rail and the middle of the clamp.

**Eyelet ring (The connection kit 094-110/CNT can be used for an easier connection to 6 mm diameter cables).

***Measured with the two clamps mounted on the guide rail.

PBH-4 & PBH-8 holders: 4 mm receptacles are available for the current (power) and voltage (sense) cables.

Prismatic and Pouch Cell Holders



PPBH-1100

4 POINT measurement

	PPBH-132	PPBH-1100
Cell min height/mm		0
Cell max height/mm		139
Min leads separation distance/mm		66
Max leads separation distance/mm		155
Number of channel		1
Max current/A	32	100
Measurement type	4 point	
Receptacles diameter/mm	4 (power and sense)	4 (power and sense) and 6 (power)
To be used with	All instruments	
Max operating T°/°C	80	
Size : HxWxD/mm	265x320x300	320x320x360
Weight/kg	3	5.1
Catalog n°	092-PC32/1	092-PC100/1

Pouch cell and prismatic holder: 4 mm receptacles can be used for currents up to 32 A. For higher currents, the 6 mm flush mounting plugs should be used. These are compatible with FlexP 0160, HCV-3048, CC4-60A and CC8 cables.



PPBH-132

4 POINT measurement

Pouch cell and prismatic holder: 4 mm receptacles are available for the current (power) and voltage (sense) cables.

Battery Accessories.

Current Collectors

Biologic's current collectors offer the possibility to connect in parallel several channels and increase the maximum current that can be passed through the cell, in order to simplify and reduce the footprint of your setup.



CC8



CC4-60A



CC4-200A

	CC4-60A	CC8	CC4-200A
Connection details			
Input			
Power cables/receptacles diameter/mm	4		6 (IP2x)
Voltage sense receptacles diameter/mm	2		4
Number of input channels	4	8	4
Max current/channel/A	15		50
Output			
Power receptacles diameter/mm	6 (IP2x)		8 (Amphenol, IP2x)
Voltage sense receptacles diameter/mm	2 (IP2x)		4 (IP2x)
Max output current/A	60	120	200
Cables details			
Output power cables	1 pair of 2 m power cables with 6 mm receptacles		1 pair of 2.5 m power cables with 8 mm receptacles and M8 threads
Output voltage cables	1 pair of 2 m sense cables with 2 mm banana plugs		1 pair of 2.5 m sense cables with 4 mm banana plugs
Instrument compatibility	BCS-815 VSP-300 VMP-300	BCS-815*	FlexP0160 FlexP0060 HCV-3048
Included connection kit	094-110/CNT**		093-200/CNT**
Measurement type	4 point		
Max operating T°/°C	80		
Size (with feet) : HxWxD/mm	70x170x88	70x300x88	120x248x169
Weight/kg	3.8		
Catalog n°	096-022	096-015/1	093-100/CC4

*The CC8 comes with BCS tablets and cables. It is also compatible with the VSP-300 and VMP-300 and can be provided without cables and tablets using the following part number: 096-015.

Sense Adapter Module (SAM-50)



SAM-50

This can be added to a multichannel system to perform stack measurements up to 60 V for 5 channel boards and a 10-element measurement. 3 SAM-50s can be linked to follow-up 30 elements.

Sense Adapter Module	Catalog n°
SAM-50	092-26

Redox Flow Battery Cells.




This range of redox flow battery cells are manufactured by **Pinflow energy storage**.

BioLogic provides two types of redox flow batteries with two different active surface areas: 5 cm² and 10 cm². This package allows you to work with both aqueous or organic electrolytes depending on your research needs or studies. For direct use, we have different testing packages available, made up of 8 felts and four membranes, that can be chosen between Fumasep F-1850, multiple Nafion types or Vanadion.


Sealings and bipolar plates can also be purchased separately as spare parts.

Complete turnkey set-ups with climatic chambers and flow control are also available. Please ask your local reseller for more information or see the table below.




The lab cells from Pinflow are specially designed to control the pressure applied on the carbon felts that are used as electrodes. Using rigid components and non-flat bipolar plates, one can not only perform reproducible experiments, but also, it is possible to use various electrode thicknesses.

The schematic below shows that by providing non-flat carbon polymer composite plates one can control the compression and the thickness of the carbon felt electrode being used. The thickness of the electrode compartment can be easily set up using carbon-polymer composite plates with a defined stump or pit.



Description		5 cm ²	10 cm ²	
Lab Cells	Standard	P-LABCELL/5	P-LABCELL/20	
	For organic media	P-LABCELL/5F	P-LABCELL/20F	
Consumables	Testing packages (8x felts, 4x membranes)	P-TEST-PACK/5-11	P-TEST-PACK/20-11	
	Nafion N115	P-TEST-PACK/5-13	P-TEST-PACK/20-13	
	Nafion N117	P-TEST-PACK/5-14	P-TEST-PACK/20-14	
	Nafion N212	P-TEST-PACK/5-15	P-TEST-PACK/20-15	
	Vanadion	P-TEST-PACK/5-16	P-TEST-PACK/20-16	
	Membranes (x1)	Fumasep F-1850	P-MEMB/5-11	P-MEMB/20-11
	Nafion N115	P-MEMB/5-13	P-MEMB/20-13	
	Nafion N117	P-MEMB/5-14	P-MEMB/20-14	
	Nafion N212	P-MEMB/5-15	P-MEMB/20-15	
	Vanadion	P-MEMB/5-16	P-MEMB/20-16	
Spare parts	Pack of sealings (x10)	Standard	P-PSEAL/5	P-PSEAL/20
		For organic media	P-PSEAL/5F	P-PSEAL/20F
	Sealings	Standard	P-SSEAL/5	P-SSEAL/20
		For organic media	P-SSEAL/5F	P-SSEAL/20F
	Bipolar plates (1x set of 2)	-	P-SBIPOL/5	P-SBIPOL/20
	Tubing	-	P-TUBINGS	
	Fittings	-	P-FPM	
Accessories	Felt cutter	-	P-CUTTER/5	P-CUTTER/20
	Wrench	-	P-WRENCH	
	Stand	-	P-STAND	
	Peristaltic pump	-	P-LABCELLPUMP/2	
	Climate chamber	-	P-CCHAMBER	
	Data acquisition module	-	P-DATAMOD	

Pinflow RFB cells exhibit low internal resistance as can be seen on the following curve. Typical values of internal resistance are lower than 2 Ω.cm².



Pinflow cells were used and characterized in the following papers:

<https://doi.org/10.1016/j.memsci.2018.02.011>

<https://doi.org/10.1016/j.jpowsour.2018.01.079>

Application note: Need more information?

In this application note, a Vanadium Redox Flow Battery (VRFB) was characterized using typical DC and AC techniques: galvanostatic charge and discharge cycling and Electrochemical Impedance Spectroscopy (EIS). [Click here](#)



Analytical Cells.

Small Volume Cells

Each voltammetry cell is designed for a specific application (specific working electrode, volume of solution, oxygen-free condition, etc.).

For example, for voltammetry investigations using standard working electrodes with an outer diameter (OD) of 6 mm, the fixed configuration of SVC-3 kit is recommended.

For applications requiring other working electrode shapes, the SVC-2 is more suitable.

If only a small amount of the electroactive compound is available, SVC-2 in the microvolume mode is recommended. Here is a list of the available cells:

- SVC-2, modular
- SVC-3, for a volume of 5 to 20 mL, only for working electrode with OD of 6 mm
- VC-4, for a volume of 1 to 3 mL, only for working electrode with OD of 6 mm
- Bulk electrolysis cell, for a volume of 100 mL



SVC-2



SVC-3






VC-4




Bulk electrolysis cell

Small Volume Cells	Catalog n°	Catalog n°	Catalog n°	Catalog n°
Products	SVC-2 A-012668	SVC-3 A-012669	VC-4 A-011224	Bulk electrolysis cell A-001197
Content				
Sample vial/mL	20 (7 pieces) A-001056	20 (7 pieces) A-001056	5 (7 pieces) A-011504	100 (1 piece) A-012632
Counter electrode (CE)/mm	57 A-002222	50 A-002222	57 A-002222	230 A-002234
PTFE cap	A-012670	A-012671	A-011226	A-012551
Purge tube (ETFE), 100 mm	-	-	-	-
Additional items	Adapter 10 to 6 mm -		Cell holder A-011227	Porous carbon electrode A-010530
				Lid for CE A-001198
				Chamber for CE A-001196
				O-ring A-001236
				Port plug A-009131
				Stirrer bar A-000178
Options				
Sample holder/mm	9.0 (2 pieces) A-012177	6.0 (2 pieces) A-012176		
Cell holder	for 20 mL A-001209	for 20 mL A-001209		
Purge tube (ETFE)/m	1 A-010537	1 A-010537	1 A-010537	1 A-010537
Working electrodes	See page 17	See page 17	See page 17	See page 17
Reference electrodes	See page 18	See page 18	See page 18	See page 18

PTFE cap




SVC-2 modularity: microvolume mode



Small amount of solution (200 µL in the sample holder)


oxygen free mode



For electrode with OD of 10 mm

For electrode with OD of 6 mm

For microelectrode with OD of 4 mm



Please note that a full, purpose-built, analytical kit is also available SK-2 (A-012763) but the reference electrode must be purchased separately.

This kit includes:

- SVC-3 kit (A-012669)
- PK-3 polishing kit (A-011975) see page 16
- one glassy carbon electrode, OD 6.0 mm, ID 3.0 mm (A-002012) see page 17
- one platinum electrode, OD 6.0 mm, ID 1.6 mm (A-002013) see page 17

Cell geometry


The geometry of the cell should be optimized to reduce the ohmic drop. Working and reference electrodes must be positioned close to one another. The counter electrode should not limit the transfer of electrons, so its contact surface should be larger than the contact surface of the working electrode.

Analytical Cells.

Large Volume Cells

This cell is perfect for standard analytical electrochemistry experiments with electrolyte volumes of several tens of mL.

It is also compatible with the BluRev product range, ie electrode rotators.




Analytical Cell kits

	EL-ELECTRO-80	EL-ELECTRO-80DJ	EL-ELECTRO-150DJ
Glass cell	EL-A-001 (80 mL)	EL-A-002 (80 mL) (double jacketed)	EL-A-020 (150 mL) (double jacketed)
PTFE cap 5 holes		EL-A-003	
PTFE ring, silicon encapsulated, OD 102 mm		EL-A-004	
Cell collar with clamp		EL-A-005	
PTFE ring, silicon encapsulated, OD 102 mm		EL-A-008	
Bridge tube for reference electrode, OD 6 mm		EL-A-009	
Purge tube	∅		EL-A-016
Double purge tube	∅		EL-A-006
Double nut 25 mm and 12 mm diameter	∅		EL-A-011
Telescopic cell stand	∅		EL-A-012
Options			
Electrode bridge extension for electroanalytical cell		EL-A-022	
Bridge tube for reference electrode of OD 8 mm		EL-A-017	
PT100 probe		EL-C-014	
220 V - Magnetic stirrer & header, without PT100		EL-C-015A	
110 V - Magnetic stirrer & header, without PT100		EL-C-015B	
Aluminum base holder for magnetic stirrer		EL-C-018	
Set of 10 porous 4 mm glass frits (CoralPor) with PTFE heat shrink (200 mm)		092-VYC4	

Multi Purpose Cells.

FlexCell®

These cells manufactured by Gaskatel are ideal for corrosion experiments in aggressive media, as well as studies on Gas Diffusion Electrodes (GDE) and membranes. Their unique design and a specific choice of materials allow repeatability and keep away from common pitfalls in other cells: heterogeneous electric field, variable ohmic drop, crevice corrosion, degradation of the cell. Used in combination with the robust HydroFlex® or MiniHydroFlex® hydrogen reference electrode, this is the perfect cell for corrosion, membrane and GDE studies. Two versions are available: PTFE and PP.



FlexCell® PTFE

FlexCell® PP

Description

Description	Catalog n°
FlexCell® PP – Electrochemical Test Cell made of PP	G-FLEXCELL/PP
FlexCell® Analyte Compartment in PP for membrane studies	G-COMP/PP
FlexCell® PTFE – Electrochemical Test Cell made of PTFE	G-FLEXCELL/PTFE
FlexCell® Analyte Compartment in PTFE for membrane studies	G-COMP/PTFE

Options to be ordered separately (see p. ???)

HydroFlex® Hydrogen Reference Electrode	G-HYDROFLEX
HydroFlex® Hydrogen Reference Electrode Starter Kit	G-HYDROFLEX-KIT
MiniHydroFlex® Hydrogen Reference Electrode	G-MINIHYDROFLEX


Operating conditions

Max. exposure time/h	24
Temperature/°C	-20 to 120°C
pH range	-2 to 16
Max. current/A	3



PTFE analyte compartment

PP analyte compartment



Cross-section of the cell and instrument connection.

Designation	Max. sample size /cm	Max. sample thickness /mm*	Active area /cm ²	Electrolyte volume /ml	Dimensions (w/o screws adapters) (H x W x D)/mm
FlexCell®	3 x 5	10	3	30	75 x 100 x 132
Analyte Compartment	3 x 5	10	3	12	75 x 100 x 30

*Longer wing screws can be supplied

Materials



Main	PP or PTFE	PP	PolyPropylene
Gasket	Silicon	PTFE	PolyTetraFluorEthylene
O rings	EPDM	PSU	PolySulfone
Cover plate	PSU	PMMA	PolyMethylMethAcrylate
Gas compartments	PMMA or PSU	EPDM	EthylenePropyleneDiene Monomer
Counter electrode	Pt-Ir	Pt	Platinum
		Ir	Iridium

Corrosion Cells.

Standard Corrosion Cells

The standard corrosion cell is available with a single or double jacket glass cell for temperature control. It is provided with two graphite rods to be used as counter electrodes, a bridge tube, to ensure the minimum distance between the reference and the working electrode, and purge tubes, to maintain a controlled gaseous environment.

Two kits are available: a standard one and an advanced kit with telescopic cell stand, sample holder and double purge tubes.




Corrosion Cell kits

	Basic corrosion cell kit EL-CORR-1	Complete corrosion cell kit EL-CORR-1DJ
Glass cell 1 L	EL-C-001	EL-C-002 (double jacketed)
PTFE cap		EL-C-003
PTFE ring, silicon encapsulated, OD 102 mm		EL-C-004
Cell collar with clamp		EL-C-005
Graphite counter electrode rod (2 pieces) $\rho = 1.070 \mu\Omega \text{ cm}$		EL-C-009
Bridge tube for 6 mm diameter reference electrode	EL-C-008	
Purge tube	EL-C-016	∅
Double purge tube	∅	EL-C-006
Double nut 25 mm and 12 mm diameter	∅	EL-C-011
Telescopic cell stand	∅	EL-C-012
Sample holder 1 cm ² (max sample thickness 3.4 mm and max diameter 14.6 mm)	∅	092-C-016
Options		
Bridge tube for 8 mm diameter reference electrode		EL-C-017
PT100 probe (indicate connector type)		EL-C-014
Magnetic stirrer & heater, without PT100 probe	220 V	EL-C-015A
	110 V	EL-C-015B
Aluminum base holder for magnetic stirrer and 1 L cell vial		EL-C-018
Set of 10 porous frits (4 mm CoralPor™) with PTFE heat shrink (200 mm)		

Avesta Cell


The Avesta Cell is an electrochemical cell developed for pitting corrosion testing (ASTM G150). It is designed to avoid microcrevice corrosion formed between the working electrode and the gasket at the bottom aperture of the cell.



Avesta Cell	Catalog n°
Avesta cell kit	EL-AVESTA
Content	
Double jacketed cell glass 250 mL	EL-AV-002
PTFE cap 5 holes	EL-AV-003
O-ring PTFE silicone encapsulated	EL-AV-004
Cell collar with clamp	EL-A-005
Double purge tube	EL-AV-006
Filter paper ring (100 pieces)	EL-AV-007
Graphite counter electrode rods (2 pieces) $\rho = 1.070 \mu\Omega \text{ cm}$	EL-C-009
Bridge tube for RE 6 mm	EL-AV-010
Double nut	EL-C-011
Telescopic cell stand	EL-AV-012
Skeleton	EL-AV-015
Options	
Peristaltic pump for low flow	EL-AV-008
Bridge tube for reference electrode with OD of 8 mm	EL-AV-013
Single purge tube	EL-AV-014
Temperature probe PT100	EL-C-014
Set of 10 porous frits (4 mm CoralPor™) with PTFE heat shrink (200 mm)	092-VYC4

A filter paper ring placed between the sample and the gasket is flooded by distilled water in order to eliminate crevice corrosion.


The water flow is controlled by a peristaltic pump (EL-AV-008).



Corrosion Cells.

Flat Cells, 1 to 10 cm² sample area

This cell with a volume of 250 mL is perfect for experiments on flat specimens of 1 or 10 cm² surface area.



This cell has a double jacket for temperature control and three holes for reference electrodes, purge tubes and temperature probes with an inner diameter of 17.6 mm and the two others of 8.3 mm.

Flat cell kit	Glass kits		PVDF kits	
	1 cm ² EL-FLAT-1CM	10 cm ² EL-FLAT-10CM	1 cm ² EL-FLAT-1CM-PVDF	10 cm ² EL-FLAT-10CM-PVDF
Platinum mesh counter electrode (54 mm wire/ 80 mesh), 25 x 35 mm	A-702439			
Reference electrode bridge tube (6 mm diameter)	EL-F-004		∅	
Glass part for flat cell (250 mL)	EL-F-002		EL-F-PVDF	
Mechanical parts	EL-FLAT-3H	EL-FLAT-4H	EL-FLAT-3H	EL-FLAT-4H

Investigations in aggressive media

If the experiment is performed in more aggressive media such as fluorhydric acid, it is possible to obtain the body of the flat cell in PVDF* instead of glass (Polyvinylidene fluoride).

	Fluorhydric acid 48%		Sulfuric acid 98%		Phosphoric acid 85%		Hydrochloric acid 35%		Nitric acid 70%		Perchloric acid		Sodium hydroxide 50%		Potassium hydroxide concentr.	
	20° C	50° C	20° C	50° C	20° C	50° C	20° C	50° C	20° C	50° C	20° C	50° C	20° C	50° C	20° C	50° C
PTFE	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
PVDF*	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Borosilicated glass	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■

Options	Catalog n°
PVDF body (single-jacketed)	EL-F-PVDF
EPDM O-rings for 1 cm ² with PEEK ferrule for CE*	EL-SEAL-1B
EPDM O-rings for 10 cm ² with PEEK ferrule for CE*	EL-SEAL-10B
PTFE O-rings for 1 cm ² with PEEK ferrule for CE*	EL-SEAL-T1B
PTFE O-rings for 10 cm ² with PEEK ferrule for CE*	EL-SEAL-T10B
EPDM O-rings for 1 cm ² *	EL-SEAL-1
EPDM O-rings for 10 cm ² *	EL-SEAL-10
PTFE O-rings for 1 cm ² *	EL-SEAL-T1
PTFE O-rings for 10 cm ² *	EL-SEAL-T10
Set of 10 porous frits (4 mm CoralPor™) with PTFE heat shrink (200 mm)	092-VYC4

* The O-ring kits include 4 O-rings for the glassware side and 10 O-rings for the sample side

Galvanic Cells, 1 to 10 cm² sample area

Thanks to the modular design of the flat cell, it is possible to place two different materials at each end of the cell.

The surface area may be 1 or 10 cm².



Sample size:
30.5 mm min. diameter (EL-GAL-1CM)
40.5 mm min. diameter (EL-GAL-10CM)
6 mm max. of thickness

Galvanic Cell kits	1 cm ² EL-GAL-1CM	10 cm ² EL-GAL-10CM
Content		
Flat cell kit 1 cm ²	EL-FLAT-1CM	EL-FLAT-10CM
Galvanic kit 1 cm ²	092-FLAT/1	092-FLAT/10

Plate Material Evaluating Cell, up to 1 cm² sample area

This cell was developed to evaluate a plate material such as metals, semi-conducting plates, etc.

The sample plate is sandwiched between the two cell blocks. The required volume of solution is about 1 mL.



Plate Material Evaluating Cell	Catalog n°
Plate material evaluating cell	A-011951
Content	
PTFE cell [body & base] (1 piece)	-
O-ring (1 piece)	-
Screw 20 mm (1 piece)	-
Purging tube, 100 mm	-
Platinum counter electrode (1 piece)	A-002222
Options	
O-ring (10 pieces)	A-012022

Coating Cell



Coating Cell	Catalog n°
Coating cell kit	EL-COAT
Content	
Glass for coating cell	EL-P-002
Nylon base with three feet	EL-P-003
Rubber cup with two holes	EL-P-004
Metallic clamp	EL-P-005
O-ring for coating cell	EL-P-006
Graphite rod counter electrode (L: 145 mm, OD: 6 mm, ρ = 1.070 μΩ cm)	EL-P-009
Options	
Bridge tube for 6 mm reference electrode	EL-P-008
Mask for 1 cm ² (20 pieces)	EL-P-011
Mask for 3 cm ² (20 pieces)	EL-P-012
Mask for 10 cm ² (20 pieces)	EL-P-013

Electrodes.

Working Electrodes

To address every application, a wide range of working electrodes (WE) is available with diameters ranging from 7 µm up to 6 mm and made of different materials.

Biologic exclusive




M-BDD-3: Boron-doped diamond 3 mm diameter disk

Standard type



Gauze type



Dimensions in mm


Working Electrodes

		Isolation	OD/mm	Electrode Ø ±4%	Catalog n°	
Carbon	Standard glassy carbon electrode	PEEK	10	5 mm	A-002417	
		PEEK	6	3 mm	A-002012	
		PEEK	6	1.6 mm	A-012297	
		PEEK	6	1 mm	A-002411	
	Small glassy carbon electrode	PEEK	3	1.6 mm	A-012298	
		PEEK	3	1 mm	A-002412	
	Standard pyrolytic graphite electrode	Basal plane	PEEK	6	3 mm	A-002252
		Edge plane	PEEK	6	3 mm	A-002253
	Standard plastic formed carbon electrode	PEEK	6	3 mm	A-002408	
		PEEK	6	1 mm	A-002409	
Small plastic formed carbon electrode	PEEK	3	1 mm	A-011854		
Platinum (99.95% purity)	Standard platinum electrode	PEEK	10	5 mm	A-002420	
		PEEK	6	3 mm	A-002422	
		PEEK	6	1.6 mm	A-002013	
Small platinum electrode	PEEK	3	1.6 mm	A-002313		
Gold	Standard gold electrode	PEEK	10	5 mm	A-002418	
		PEEK	6	3 mm	A-002421	
		PEEK	6	1.6 mm	A-002014	
Small gold electrode	PEEK	3	1.6 mm	A-002314		
Silver	Standard silver electrode	PEEK	10	5 mm	A-002416	
		PEEK	6	3 mm	A-002419	
		PEEK	6	1.6 mm	A-002011	
Small silver electrode	PEEK	3	1.6 mm	A-002315		
Palladium	Standard palladium electrode	PEEK	6	1.6 mm	A-002019	
	Small palladium electrode	PEEK	3	1.6 mm	A-002319	
Nickel	Standard nickel electrode	PEEK	6	1.5 mm	A-002016	
Copper	Standard copper electrode	PEEK	6	1.6 mm	A-002017	
		PEEK	6	3 mm	A-012584	
Iron (99.65% purity)	Standard iron electrode	PEEK	6	1.5 mm	A-002018	
		PEEK	6	3 mm	A-012585	
Carbon paste	Standard carbon paste electrode hole depth 4 mm	PEEK	6	3 mm	A-002210	
	Small carbon paste electrode hole depth 4 mm	PEEK	3	1.6 mm	A-002223	
	Carbon paste oil base 1 g				A-001010	
Boron-doped diamond	Doping level between 500 and 1000 ppm. The electrode is a 500 µm thick disk attached to a conductive rod in brass. It is polished with an Ra < 10 nm.	PEEK	7	3 mm	M-BDD-3	


Maintenance of working electrodes

To refresh the electrode surface, we recommend polishing it before each measurement.

PK-3 electrode polishing kit	Catalog n°
Content	
0.05 µm polishing alumina (20 mL)	A-001050
1 µm polishing diamond (10 mL)	A-002054
Glass plate (1 piece)	A-002249
Alumina polishing pad (10 pieces)	-
Diamond polishing pad (10 pieces)	-
Spare parts	
Alumina polishing pad (20 pieces)	A-001040
Diamond polishing pad (20 pieces)	A-001041
Emery paper UF800 (20 pieces)	A-012611
Coarse polishing pad (20 pieces)	A-001042
6 µm polishing diamond (10 mL)	A-002053



Counter Electrodes




Counter Electrodes	Size/mm	Wire Ø/mm	Surf. area/cm²	Purpose	Catalog n°
Platinum*	57	0.5	~ 0.7	SVC-2, VC-4, plate material evaluating cell	A-002222
	50	0.5	~ 0.7	SVC-3	A-002233
	230	0.5	~ 3.6	RRDE-3A, bulk electrolysis, SVC-3	A-002234
Gold	230	0.5	~ 3.6	RRDE-3A, bulk electrolysis, SVC-3	A-012638
Nickel	230	0.5	~ 3.6	RRDE-3A, bulk electrolysis, SVC-3	A-012639
Platinum* gauze electrode, PEEK body	25x35	0.08	~ 22.9		A-002250
Platinum* gauze electrode, 54 mm wire	25x35	0.08	~ 22.9	Flat cell	A-702439
Platinum* electrode for bulk electrolysis cell	40x50	0.5	~ 47.4	Bulk electrolysis cell	A-013612
Gold gauze electrode, PEEK body	25x35	0.07	~ 29		A-002251

Electrodes.

Small-Size Reference Electrodes for aqueous media

Reference electrodes are divided into two groups according to the media in which the electrode is immersed (aqueous or organic media).



Small-Size Reference Electrode for Aqueous Media (Ag/AgCl, Hg)

	Junction	Electrolyte	Purpose	Catalog n°
RE-1B Ag/AgCl reference electrode ⁽¹⁾	IPPG*	3 M NaCl	SVC-2, SVC-3, VC-4, bulk electrolysis, RRDE-3A, flat cell	A-012167
RE-1S Ag/AgCl reference electrode ⁽¹⁾	IPPG*	3 M NaCl	SECM	A-012168
RE-1BP reference electrode (Ag/AgCl)	Ceramic	3 M NaCl	SVC-2, SVC-3, VC-4, bulk electrolysis, RRDE-3A, EQCM, flat cell	A-013613
RE-1CP Ag/AgCl reference electrode	Ceramic	Saturated KCl	SVC-2, SVC-3, VC-4, bulk electrolysis, RRDE-3A, flat cell	A-013429
RE-2BP Hg/Hg ₂ Cl ₂ reference electrode	Ceramic	Saturated KCl	SVC-2, SVC-3, VC-4, bulk electrolysis, RRDE-3A, flat cell	A-013430
RE-2CP Hg/Hg ₂ SO ₄ reference electrode, free from chloride	Ceramic	Saturated K ₂ SO ₄	SVC-2, SVC-3, VC-4, bulk electrolysis, RRDE-3A, flat cell	A-013431
RE-61AP Hg/HgO reference electrode main body in polyacetal resin	Ceramic	1 M NaOH	For alkaline media	A-013395

Spare parts


Set of 10 porous 2.8 mm glass frits (CoralPor™) with high temperature PTFE heat shrink (200 mm long, 3.2 mm diameter). Only compatible with: A-012167, A-01277, A-012178	092-VYC3
Set of 10 porous 2.8 mm glass frits (CoralPor™) with low temperature PVDF heat shrink (200 mm long, 3.2 mm diameter). Only compatible with: A-012167, A-01277, A-012178	092-VYC3A
Set of 10 porous 2.8 mm glass frits (CoralPor™) with PTFE heat shrink (200 mm long, 4.8 mm diameter) Only compatible with A-012168	092-VYC5

Options

RE-PV preservative vial for reference electrode, 10 mL	A-012108
Bridge tube Ø 9.0 mm (2 pieces)	A-012177
Bridge tube Ø 9.0 mm (22 pieces)	A-012307

* Ion Permeable Porous Glass

Small-Size Reference Electrodes for non aqueous media



Small-Size Reference Electrodes for Non Aqueous Media (Ag/Ag⁺)

	Junction	Electrolyte	Purpose	Catalog n°
RE-7 non aqueous reference electrode (Ag/Ag ⁺)	IPPG*	Ag/Ag ⁺ /ACN**/ TBAP***	CV	A-012171
RE-7S non aqueous reference electrode (Ag/Ag ⁺)	IPPG*	Ag/Ag ⁺ /ACN**/ TBAP***	SECM	A-012172

Spare parts

Electrolyte solution (10 mL)	A-012549
PTFE cap with Ag wire (for RE-7)	A-012057
Sample holder 6 mm diameter (for RE-7) (2 pieces)	A-012176
Set of 10 porous 2.8 mm glass frits (CoralPor) with high temperature PTFE heat shrink (200 mm long, 3.2 mm diameter). Only compatible with: A-012171, A-012177, A-012178	092-VYC3
Set of 10 porous 2.8 mm glass frits (CoralPor) with low temperature PVDF heat shrink (200 mm long, 3.2 mm diameter). Only compatible with: A-012171, A-012177, A-012178	092-VYC3A
Set of 10 porous 2.8 mm glass frits (CoralPor™) with PTFE heat shrink (200 mm long, 4.8 mm diameter) Only compatible with A-012171	092-VYC5

Options

RE-PV preservative vial for reference electrode, 10 mL	A-012108
Bridge tube Ø 9.0 mm (2 pieces)	A-012177
Bridge tube Ø 9.0 mm (22 pieces)	A-012307

- * Ion Permeable Porous Glass
- ** Acetonitrile
- *** Tetra Butyl Ammonium Perchlorate

Don't forget!

There is a huge amount of supporting information on www.biologic.net. The BioLogic Learning Center has over 130 articles and if you need to drill down deeper, there are over 80 application notes and 50 technical notes for the field of electrochemistry alone.

www.biologic.net

Support: How to check your reference electrode

i BioLogic's Learning Center is a great source of information for tech-tips, theory and product information. Visit our Learning Center or [click here](#) to find out how why reference electrode maintenance is so important.

Maintenance of Reference Electrodes.

Store your reference electrode immersed in the electrolyte

When not in use, we recommend that you keep reference electrodes in sealed, air-tight vials in order to prolong their life. The storage solution should be identical to the filling solution of the reference electrode. Prefer a cold and dark place.

Prevent contamination

To prevent contamination of the reference electrode, a sample holder can be used.

Replace the junction when needed

If you are using IPPG junctions, yellowish discoloration indicates contamination. This is caused by the absorption of organic compounds into the pores. The average pore diameter of IPPG is about 40 - 200 Å. If you are using CoralPor™ junctions, you might want to use one of the available replacement kits. The average pore diameter of CoralPor™ is about 4 - 10 nm.

Potentials of common reference electrodes

E/V vs. NHE at 25 °C	
0.930	Hg/HgO/ NaOH (0.1 M)
0.650	Hg/Hg ₂ SO ₄ / K ₂ SO ₄ (sat)
0.624	Fc/Fc ⁺ TBAP (0.1M) ACN
0.542	Ag/Ag ⁺ TBAP (0.1M) ACN
0.241	Hg/Hg ₂ Cl ₂ KCl (sat)*
0.236	Hg/Hg ₂ Cl ₂ NaCl (sat)*
0.205	Ag/AgCl/ KCl (3.5 M)
0.197	Ag/AgCl/ KCl (sat)
0.194	Ag/AgCl/ NaCl (sat)
0.000	NHE Normal Hydrogen Electrode

* Hg/Hg₂Cl₂ : Calomel

Hydrogen Reference Electrodes.


These hydrogen reference electrodes, manufactured by Gaskatel, are beneficial in that they are Hg-free and not made of glass, which extends their range of operating conditions. These robust reference electrodes are easy to use and robust. Hydrogen is contained within a cartridge that is easily replaceable.

Description	Catalog n°
HydroFlex® Hydrogen Reference Electrode	G-HYDROFLEX
HydroFlex® Starter Kit*	G-HYDROFLEX-KIT
MiniHydroFlex® Hydrogen Reference Electrode	G-HYDROFLEX-CARTRIDGE

Options	
HydroFlex® Hydrogen Cartridge (x4)	G-HYDROFLEX-CARTRIDGE
MiniHydroFlex® Hydrogen Cell (x2)	G-MINIHYDRO-CELL
Bridge tube for use of HydroFlex® with EL-ELECTRO	EL-A-017
Bridge tube for use of HydroFlex® with EL-CORR	EL-C-017
Bridge tube for use of HydroFlex® with EL-FLAT	EL-F-004B

Operating conditions	
Cartridge lifetime/months	12 (MiniHydroFlex®), 6 (HydroFlex®)
Temperature/°C	-20 to 120 (PTFE body)
pH range	-2 to 16

The HydroFlex® electrodes are compatible with the FlexCell®.



Replacement kits	Compatible reference electrodes	Compatible bridge tubes	Content
092-VYC3	A-012167 A-012178	A-012176 A-012306 A-012177 A-012307	10 glass frits (φ 2.8mm Coralpor™) 200 mm long heat shrink tube (φ 3.2 mm)
092-VYC4	∅	EL-C-005 EL-C-017 EL-F-004B EL-F-004 EL-A-017 EL-A-008	10 glass frits (φ 4 mm CoralPor™) 200mm long heat shrink tube (φ 4.8 mm)
092-VYC5	A-012168	∅	10 glass frits (φ 2.8 mm CoralPor™) 200 mm long heat shrink tube (φ 4.8 mm)



Part	Materials
Body	PTFE, PP
Shaft	PEEK, PTFE
Cap and cartridge	PC or PVC
Measuring electrode	Pt, Pd
PP	PolyPropylene
PTFE	PolyTetraFluorEthylene
PVC	PolyVinylChloride
PC	PolyCarbonate
PEEK	PolyEtherEtherKetone
Pt	Platinum
Pd	Palladium

HydroFlex® Starter Kit Content	
1 HydroFlex® Hydrogen reference electrode incl. 1 hydrogen cartridge	
1 connector lead	
1 hydrogen cartridge	
1 cartridge wrench	
1 operating time wrench	

Connection Accessories.

High-Temperature Extension Cables

These extension cables are intended to be connected between the instrument's cell cable and the cell.

Available in two different sizes, 1.3 m and 2.5 m long and with a maximum temperature of 150 °C, these cables are compatible with:

SP-50e, SP-150e, VSP, VSP-3e, VMP-3e and boosters for Essential range instruments (20 A maximum)
 SP-200, SP-240, SP-300, VSP-300, VMP-300, and boosters for Premium range instruments (20 A maximum) BCS-805, BCS-810 and BCS-815

High-temperature cables	1.3 m	2.5 m m
Temperature /°C	-40 °C to 150 °C	
Max current /A (2 mm connectors)	2	
Max current /A (4 mm connectors)	20	
Cable diameter /mm	12.7	
Catalog n°	092-25/1	092-25/2



Glove Box Cables



Hermetic cell cable for glove box

As standard, the potentiostat and the booster are provided with a 1.5 m long cell cable. The cable connected from the booster to the potentiostat is 0.8 m long for VMP3 based instruments.

For some applications, the user may need different length cables. For this reason, longer cables are available (for more information, contact your sales representative).

For applications carried out in glove boxes, cell cables are also available.

Hermetic cell cable for glove box		
	VMP3	VMP-300
Catalog n°	092-23/5	094-101/6 (standard cable) 094-101/8 (low current cable)
Content:		
Feedthrough type	12 pins	25 pins*
Inside glovebox (length: 1 m)	Cable with 2 mm connectors on one side and 12-pin Jaeger connector on the other side	Cable with electrometer on one side and 25-pin connector on the other side
Outside glovebox (length: 1.5 m)	Cable with SubD25 connector on one side and 12 -pin Jaeger connector on the other side (length 1.5 m)	Cable with SubD25 connector on one side and 12 -pin Jaeger connector on the other side (length 1.5 m)
Requirement:		
Hole to make in the glove box/mm	27	45

*Two feedthrough seals one installed in the glove box wall the other dedicated to the channel board calibration outside the box

Set-up connection

Bad connections can affect measurements (stability of potentiostat, artefacts etc).

In order to optimise your set-up, we recommend you use the accessories described in this section.

Multi-electrode investigation cables

For the Essential product range, we offer several options to facilitate the use of connection cables when multi-electrode experiments are performed, for example, RRDE experiments or corrosion experiments on several samples using the same reference and the same counter electrode.



Nstat box (8 channels)



Bipot cable: dedicated to RRDE applications

Description	Nb of channels	Length	Catalog n°
Bipot cable (for SP-150e, VSP, VSP-3e and VMP-3e)	2	1.5 m	092-22/12
Nstat box (for VSP, VSP-3e and VMP-3e). External power supply required.	4	1.5 m	092-16
	4		092-22/3
External power supply for the Nstat box			096-16/1

Connectors

	Content	Catalog n°
Alligator clips	2 mm 8x Red, 8x Blue, 8x White, 4x Black	092-1001/40
	4 mm 3x Red, 3x Blue, 3x Black	092-1001/41
Receptacles	2 mm 20x Red, 20x Blue, 20x White, 15x Black	092-1001/42
	4 mm 12x Red, 12x Blue, 12x Black	092-1001/43
Banana plugs	2 mm 10x Red, 10x Blue, 10x White, 10x Black	092-1001/44
	4 mm 8x Red, 8x Blue, 8x White	092-1001/45
Adaptors	2mm recep.-4mm banana plugs 16x Red, 16x Blue	092-1001/47
	4mm recep. to 2mm banana plugs (x20)	092-1001/48

Connection kits

Connection kits	Catalog n°
For standard board 4 alligator clips of 2 mm: red, blue, white, black 3 receptacles of 2 mm: red, blue, white	092-1001/30
For booster board 3 alligator clips of 2 mm: red, blue, white 2 alligator clips of 4 mm: red, black 3 receptacles of 2 mm: red, blue, white 2 receptacles of 4 mm: blue, white	092-1001/31
For HCV-3048, FlexP 0060, 0160, CC4-60A and CC8 Contains: 2x6 mm receptacles: red, blue 4 lugs with 4 mm receptacles: red, blue	094-110/CNT
For FlexP0012 and CC4-200A Contains: 2 Amphenol 8 mm receptacles, red, black 4 lugs with 4 mm receptacles red, blue	093-200/CNT



094-110/CNT



093-200/CNT

Connection Accessories.

External Device Connection



IS1



DB9-8 BNC

External device connection		Catalog n°
DB9-8 BNC connector for auxiliary I/O		092-22/1
IS1 isolation module for auxiliary I/O for VMP-300 based instruments		094-081/5
Pt Probe		Catalog n°
PT100 temperature probe, to be connected to the auxiliary I/O, temperature range: -50 °C to 250 °C, Dimensions: 3 x 20 mm, Length of cable: 2.5 m, Accuracy: ±1 °C		092-22/13
For T° measurement in air		
PT100 temperature probe for T° measurement in solution with SubD9 connector		EL-C-014
PT100 temperature probe for T° measurement in solution with triad connector		EL-C-014/1

Test Boxes

Test Boxes		Catalog n°
Test Box 2	Several circuits with high precision resistors, for calibration and validation	092-22/6
Test Box 3	Three circuits: linear, two non-linear systems (Tafel & passivating) for teaching and demonstration	092-22/7



Test Box 2



Test Box 3

Faraday Cages


To avoid any external perturbations, especially for low current applications, we recommend using a Faraday cage.

Please note that in order to ensure that the cage is fully functional, it must be earthed by connecting it to the ground (this is done via a green plug on the instrument's rear panel).



Faraday cages		Catalog n°
FC-45 Faraday cage, 450x450x450 mm		094-084/1
Stand for FC-45		094-084/2

Spectroelectrochemistry.



Spectrometer System


Spectroelectrochemistry (SEC) can be useful to elucidate electrochemical reaction mechanisms. The spectroelectrochemical kit is made up of three parts (spectrometer, light source and cuvette holder).

The spectrometer is equipped with a trigger to synchronize electrochemical and spectroscopic measurements.

Spectrometer specifications


Detector	2048 element linear silicon CCD array
Full description	SEC2021-025-DUVN
Detector range/nm	200 - 1025
Grating	Blaze wavelength (300 nm)
Slit/μm	25
Wavelength resolution/nm	1.3
Fiber connector	SMA905 Core diameter: 600 μm NA=0.22
Interface	USB2.0
Operating system	Windows™ 7 / 8.1/10 (32bit / 64bit)
Dimensions (HxWxD)/mm	32x86x110

SEC2020 Spectrometer system		Catalog n°
SEC2020 spectrometer kit		A-013609
Content		
SEC2021 Spectrometer (x 1)		-
SEC2022 Deuterium halogen light source (x 1)		-
SEC2023 Cuvette holder (x 1)		-
SEC2024 Platform (x 1)		-
AC adaptor (x 1)		-
Power cable (x 1)		-
USB cable (x 1)		-
Collimator (x 2)		-
Fiber collimator (x 1)		-
Platform screw (x 7)		-
External device connection trigger cable (x 1)		-
Light source control trigger cable (x 1)		-
Plastic cuvette (x 1)		-
SMA905 adaptor for light source (x 3)		-
SMA905 adaptor for light shielding (x 2)		-
Hexagon wrench 0.89 mm (x 1)		-
Hexagon wrench 1.50 mm (x 1)		-
Software (USB memory) (x 1)		-
Waterproof box (x 1)		-
Quick manual, wavelength calibration data sheet, linearity test data sheet and warranty certificate are also included.		
Option		
Connecting cable to synchronize the SEC2020 with BioLogic instrument.		092-22/11



The SEC2020 spectrometer system uses the Czerny-Turner optical mount. This system is an M-shaped structure symmetrical to the grating (4) and is an optical system with extremely small aberration.

Light source structure




Light source specifications	
Light type	Deuterium halogen light source
Wavelength range/nm	200 - 1700
Stability	<0.1%
Drift/h	0.25%
Bulb life/h	>1000 (D2 lamp) >2000 (halogen lamp)
Fiber connector	SMA905
Size (HxWxD)/mm	46x100x165

- 1. SMA905 Connector
- 2. Slit
- 3. Collimating mirror
- 4. Grating
- 5. Focus mirror
- 6. 2048 element CCD array




Utilization modes

Transmittance	
Absorbance/transmittance	- Concentration of chemicals (solution) - Polymer extrusion processes - DNA quantification
Reflectance	- Freshness testing - Film thickness/composition (quality control) - Activation energy of photocatalytic species - Textile quality control
Fluorescence	
Fluorescence	- Marine organisms - Biology (DNA, protein, cell proliferation assay, histamine-analysis, alga monitoring) - Environmental fields (waste water analysis, ground water trace studies, hydrocarbon detection, dissolved oxygen) - Plant efficiency (plant physiology, plant breeding, horticulture, agronomy, agrochemicals, pollution studies, forestry, ecology) - Tissue diagnosis
Scattering	- Oil concentrations of oil/water system - Raman spectroscopy - Physical transition phenomena (e.g. melting point, glass transition crystallize temperature)
Irradiance	
Emission	- Astronomy (e.g. spectra of Hale-Bopp, plasma monitoring) - In situ metal monitoring - Luminescence (PL,EL), LED & laser wavelength

Static Cell



Quartz glass spectroelectrochemical cell kit	Platinum	Catalog n°	Gold	Catalog n°
0.5 mm		A-012813		A-012814
Content				
Platinum counter electrode		A-012609		A-012609
Thin layer quartz glass cell		A-012815		A-012815
PTFE cap		A-011501		A-011501
Purging tube (ETFE, 100 mm)		-		-
Gauze working electrode	80 mesh, height 5 mm	A-012606	100 mesh, height 5 mm	A-012607
1 mm		A-013510		A-013511
Content				
Platinum counter electrode		A-012906		A-012906
Thin layer quartz glass cell		A-012907		A-012907
PTFE cap		A-011501		A-011501
Purging tube (ETFE, 100 mm)		-		-
Gauze working electrode	80 mesh, height 6 mm	A-011498	100 mesh, height 6 mm	A-012017
Options				
RE-1BP Ag/AgCl reference electrode				A-013613
RE-7 non aqueous reference electrode				A-012171
Purging tube (ETFE), 1 m				A-010537



Scanning Product Accessories.

Probes

A range of probes dedicated for use with our SECM, SVP, SKP and LEIS scanning probe applications are available for the M370 and M470 systems. SECM probes can also be used with SECM150



Probes	Catalog n°
Fused silica-based SECM 10 µm diameter Platinum disk	U-23/10
Fused silica-based SECM 15 µm diameter Platinum disk	U-23/15
Fused silica-based SECM 25 µm diameter Platinum disk	U-23/25
Capillary based SECM 1 µm diameter Platinum disk	U-P5/1 ^{1 2}
Capillary based SECM 2 µm diameter Platinum disk	U-P5/2 ^{1 2}
Capillary based SECM 5 µm diameter Platinum disk	U-P5/5 ¹
Capillary based SECM 10 µm diameter Platinum disk	U-P5/10
Capillary based SECM 15 µm diameter Platinum disk	U-P5/15
Capillary based SECM 25 µm diameter Platinum disk	U-P5/25
SKP 500 µm diameter	U-SKP370/1
HR SKP 150 µm diameter	U-SKP-150
LEIS	U-LEIS370/1
SVP	U-SVP370/1

¹Unsuitable for ic-SECM/ ²Unsuitable for M370

VCAM3 Video Microscope System

The VCAM3 is a long working distance video microscope which allows users to view the positioning between probe tip and sample surface in many scanning probe electrochemistry techniques.




VCAM3 specifications

Working Distance	130 mm
Min illumination/lux	0.003
Field of view/mm	1.6*2 mm (max) to 13*8mm (min)
Operation temperature/°C	-30 to +70
Catalog n°	U-VCAM3



M470

Local Hero



If you need a local view of electrochemistry, look no further than the BioLogic M470 scanning electrochemical workstation.

A fully modular instrument, the M470 can be tailored to scan 9 techniques from 7 modules. And user-friendly software means that no time will be lost learning to operate the instrument, whether you are a specialist in local electrochemistry, or new to the field.

Techniques Available

SECM	SDS
ic-SECM	ac-SDS
ac-SECM	SKP
LEIS	OSP
SVP	



Visit www.biologic.net/products/m470/ or click [here](#) for more information

Scanning Product Accessories.

Cells

Three cells are available :
 The TriCell™ is a large volume, wide scan range cell, dedicated to LEIS, SVP, SKP, SDS techniques.
 The μTriCell™ and its Shallow version are dedicated to SECM techniques (dc, ac and ic mode). The Shallow μTriCell™ contains a slightly smaller volume of electrolyte than the μTriCell™ is more accessible and ideal for ic-SECM.


The Foil Cell has been designed for use with flat, foil type samples, such as those used for battery electrodes. It has been designed to mount directly on the baseplate of the μTriCell™ and Shallow μTriCell™.

All TriCells accommodate samples mounted in a 32 mm diameter resin cylinder.

Cells	Volume (mL)	Catalog n°
TriCell™	700	U-TRICELL
μTriCell™	7	U-uTRICELL
Shallow μTriCell™	6	U-SuTRICELL
Foil Cell	1	U-uFoilCell




TriCell™



μTriCell™



Shallow μTriCell™



Foil Cell™

M470 Glovebox Cables

The M470 Scanning Electrochemical Workstation is supplied with a full set of standard cables for use with all techniques.

For applications requiring use of the M470 in a glove box, additional cell cables are available.

These are supplied as a set of internal, feed-through and external cables to replace a single cable.



Hermetic scan stage cable for glove box

	Electrometer	Piezo Strain Gauge	Piezo Drive	Scan Stage	3300
Content :					
Feedthrough Type/pins	8	6	3	8	25
Inside glove box	Electrometer cable connects directly to 8 pin LEMO feedthrough.	Cable with connector to piezo strain gauge on one side and 6-pin Jaeger connector on the other side (length 1.5 m)	Cable with connector to piezo drive on one side and 3-pin Jaeger connector on the other side (length 1.5 m)	Cable with connector to scan stage on one side and 8-pin Jaeger connector on the other side (length 1.1 m)	Cable with 4 mm connectors on one side and 25-pin Jaeger connector on the other side (length 1 m)
Outside glove box	1 m cable with connector to SCV470 in one side and 8 pins LEMO connector in the other side	Cable with connector to SCV470 on one side and 6-pin Jaeger connector on the other side (length 1 m)	Cable with connector to SCV470 on one side and 3-pin Jaeger connector on the other side (length 1 m)	Cable with connector to SCV470 on one side and 8-pin Jaeger connector on the other side (length 1.1 m)	Cable with connector to SCV470 on one side and 25-pin Jaeger connector on the other side (length 1 m)
Requirement:					
Hole to make in the glove box/mm	12.1	21	21	27	45
Max. Required	1	1	1	3	2
Catalog n°	U-HC470ELE-L	U-HC470PSG	U-HC470PD	U-HC470STG	U-HC3300CL

*When used with an M470 and SP-300 configuration the SP-300 Ultra Low Current (ULC) glove box cable sets are also required.

USB-PIO

The USB-PIO, designed for use with the M470, allows external devices to be switched on and read. It can be used to control up to four different channels individually or collectively using the M470 software.

The USB-PIO can interface directly to user supplied cables, or to the supplied breakout PCB using the DB25 pin female connector.




USB-PIO Specifications

Temperature/°C	0 - 50
Power supply	100 VAC-240 VAC, 50/60 Hz
Fuse/A	1
Gas pressure	< 34 kPa max.
Interface	IEEE-488, RS-232C
Weight/kg	3.8

Material Testing Accessories.

How to identify the correct solution

As the electrical properties of materials depend on temperature, a Temperature Control Unit and a Sample Holder are needed to control the temperature and hold the sample (solid, pasty and liquid sample) between the parallel plates.



Temperature Control Units

High Temperature Furnace (HTF-1100)

HTF-1100 is a horizontal laboratory tube furnace dedicated to the electrical characterization of materials and to heat treatment in the temperature range between the ambient and 1100 °C.



HTF-1100 & HTSH-1100

Intermediate Temperature System (ITS-e)

ITS-e is a compact temperature chamber dedicated to the characterization of electrical properties of materials by impedance spectroscopy under controlled atmosphere in the temperature range between -35 °C and 150 °C.



ITS

Temperature Control Unit	Operating Temp.	Features	Catalog n°
HTF-1100	RT to 1100 °C	Heating rate adjustable K-type thermocouple	097-110
In-plane ITS-e	-35 to 150 °C	Temperature accuracy: 0.3 °C	097-140e
Through-plane ITS-e	-35 to 150 °C	PT1000 probes	097-140e

Sample Holders

High Temperature Sample Holder (HTSH-1100)

HTSH-1100 is a test fixture dedicated to the characterization of electrical properties of materials at high temperatures (up to 1100°C).



HTSH-1100

Controlled Environment Sample Holder (CESH-e)


CESH-e is a sample holder designed to the electrical characterization of flat material samples in parallel plate configuration in the temperature range between -40 °C and 150 °C.



CESH-e on its base

High Temperature Conductivity Cell (HTCC)

HTCC is a parallel plate conductivity cell dedicated to laboratory measurements in a wide conductivity range (200 mS/cm to 2 μS/cm).



HTCC

Other accessories	Catalog n°
Cable kit for CESH-e and MTZ-35	097-150e/01
Cable kit for CESH-e and Potentiostat	097-150e/02
Thickness measurement kit	097-150e/10
Liquid conductivity kit	097-160

Sample Holders	Operating Temp.	Features	Compatibility	Catalog n°
HTSH-1100	RT to 1100°C	Φ=25 mm	Quartz tube for controlled atmosphere Leak-tight up to 2 bar relative K-type thermocouple	HTF-1100
		Φ=12 mm		Tubular furnaces
		Φ=03 mm		
CESH-e	-40 to 150 °C	In-plane	Leak-tight up to 2 bar relative	ITS
		Through-plane		Other temperature units
HTCC	-50 to 180 °C	Platinized (x1)	Cell factor: K = 1 +/- 5% cm ⁻¹ volume: 0.5 - 1.0 mL	MCS 10
		Non-platinized (x1)		



Premium.

When only the **best** will do

Essential.

Tools for **Electrochemists**



VMP-300 16 channels,
150 A down to 100 fA



VMP-3e, 16 channels,
± 1 A down to 20 nA



SP-200
1 Channel, 500 mA
down to 100 fA



SP-240
1 Channel, ± 4 A
down to 100 fA



SP-300
2 Channels, 10 A
down to 100 fA



VSP-300
6 Channels, 40 A
down to 100 fA



SP-50e
1 channel, ± 1 A
down to 20 nA



SP-150e
2 channels, 1 A
down to 20 nA



VSP-3e
8 channels,
± 1 A down to 20 nA

BioLogic Premium is a range of state-of-the-art potentiostats designed for researchers who need the fastest, most precise, potentiostats available.

Premium potentiostats boast some of the most powerful specifications available: 100 fA to 150 Amps, 7 MHz EIS measurements and a sampling rate that can reach 1 data point every μ s.

Premium by name. Premium by nature.

[Click here for more details](#)

Shaping the future. Together.

A range of powerful, modular, high-precision, research-grade potentiostats built to handle almost any academic or industrial application imaginable.

From 1 to 16 channels. 1 Amp (native), up to 800A with boosters, EIS and Quality Indicators.


Essential measurement tools. Whatever your area of specialization.

[Click here for more details](#)

www.biologic.net

BioLogic Science Instruments

serves its customers worldwide through its affiliate offices and an extensive distribution network



Headquarters

BioLogic SAS

4, rue de Vaucanson - 38 170 Seyssinet-Pariset
France

Phone: +33 476 98 68 31

Fax: +33 476 98 69 09

Affiliate offices

BioLogic USA, LLC

P.O.Box 30009 - 37 930 Knoxville
USA

Phone: +1 865 769 3800

Fax: +1 865 769 3801

BioLogic Pvt Ltd

Unit No.204, Odyssey IT Park, Road No. 9, MIDC
Wagle Estate, Thane West, Mumbai-400604 MH
India

Phone: +91 022 25829074

BioLogic Middle East

DMCC Office - No. 1206, Indigo Icon Tower, Cluster F,
Dubai JLT

United Arab Emirates

Phone: +971 455 151 82

Fax: +971 455 151 63

BioLogic Singapore

No.18, Boon Lay Way #02-105 (A), Tradehub 21,
Singapore

Phone: +65 92335838

BioLogic Science Instruments Ltd

Bradbury House, Market Street, Glossop
Derbyshire, SK13 8AR

United Kingdom

Phone: +44 (0) 333 012 4056

BioLogic Spain

Calle Diputació, 262, 11, Barcelona 08007
Spain

Phone: +34 656 876773

BioLogic Science Instruments GmbH

Rodeweg 20, D-37081 Göttingen
Germany

Phone: +49 (0)551 38266900